




# Recombinant Human Secretogranin-1 (CHGB), partial

<b>Product Code</b>	CSB-EP005345HU1
<b>Abbreviation</b>	Recombinant Human CHGB protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P05060
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	HSTHYRASEEEPEYGEEIKGYPGVQAPEDLEWERYRGRGSEEYRAPRPQSE ESWDEEDKRNYPSELDKMAHGYGEESEEEERGLEPGKGRHHRGRGGEPRA YFMSDTREEKRFLGEGHHRVQENQMDKARRHPQGAWKELDRNYLNYGEEG APGKWQQQGDLDQTKENREEARFQDKQYSSHHTAEKRKRLGELFNPYYDPL QWKSSHFERDNDNFLEGEENELTLNEKNFFPEYNYDWWEEKPFSEDV NWGYEKRNLARVPKLDLKRQYDRVAQLDQLLHYRKKSAEFPDFYDSEEPVST HQAENEKDRADQTVLTEDEKKELENLAAMDLELQKIAEKFSQRG
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	CHGB
<b>Expression Region</b>	328-677aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-GST-tagged
<b>Mol. Weight</b>	73.5 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Description

Recombinant Human CHGB protein expression is initiated by transfecting e.coli cells with a DNA expression vector containing the gene encoding the protein of interest (328-677aa). The subsequent culturing of the cells induces the expression of the desired protein. The protein is equipped with a N-terminal 6xHis-GST tag. The resulting recombinant Human CHGB protein is collected and purified from the cell lysate through affinity purification, demonstrating a purity greater than 85%, as indicated by SDS-PAGE.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.