




# Recombinant Human Cartilage intermediate layer protein 2 (CILP2), partial

<b>Product Code</b>	CSB-EP005438HU
<b>Abbreviation</b>	Recombinant Human CILP2 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IUL8
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	HNAGGSHPRTRGQLYGLRDARSVRDPERPGTSAACVEFKCSGMLFDQRQVD RTLVTIMPQGSCRRVAVNGLLRDYLTRHPPPVAEDPAAFSMLAPLDPLGHNY GVYTVTDQSPRLAKEIAIGRCFDGSSDGFSREMKADAGTAVTFQCREPPAGR PS
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	CILP2
<b>Expression Region</b>	956-1113AA
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-GST-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	47.6 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Description

In the process of producing the recombinant Human CILP2 protein, e.coli cells are transfected with a DNA expression vector containing the gene for the desired protein (956-1113AA). Subsequent culturing induces the expression of the intended protein. The protein is fused with a N-terminal 10xHis-GST tag and C-terminal Myc tag. The recombinant Human CILP2 protein is then collected and purified from the cell lysate through affinity purification, displaying a purity exceeding 85%, as verified by SDS-PAGE.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.