




Recombinant Human Clathrin heavy chain 2 (CLTCL1), partial

Product Code	CSB-EP005594HU
Relevance	Clathrin is the major protein of the polyhedral coat of coated pits and vesicles. Two different adapter protein complexes link the clathrin lattice either to the plasma mbrane or to the trans-Golgi network .
Abbreviation	Recombinant Human CLTCL1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53675
Alias	Clathrin heavy chain on chromosome 22 ;CLH-22
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	LLVLSPRLDHTWTVSFFSKAGQLPLVKPYLRSVQSHNNKSVNEALNHLLTEEE DYQGLRASIDAYDNFDNISLAQQLEKHQLMEFRCIAAYLYKGNWWAQSVEL CKKDHLKYDAMQHAESRDAELAQKLLQWFLEEGKRECF
Research Area	Signal Transduction
Source	E.coli
Target Names	CLTCL1
Expression Region	1423-1566aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.