




Recombinant Human Cytochrome P450 2J2 (CYP2J2)

Product Code	CSB-EP006430HU
Abbreviation	Recombinant Human CYP2J2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51589
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Cytochrome P450 2J2 (CYP2J2)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	DFLKRRRPKNYPPGPWRLPFLGNFFLVDFEQSHLEVQLFVKKYGNLFSLELGD ISAVLITGLPLIKEALIHMDQNFGNRPVTPMREHIFKKNGLIMSSGQAWKEQRR FTLTALRNFGGLGKKSLEERIQEEAQHLTEAIKEENGQPFDPHFKINNAVSNII CSI TFGERFEYQDSWFQQLLKLLDEVTYLEASKTCQLYNVFPWIMKFLPGPHQTLF SNWKKLKL FVSHMIDKHRKDOWNPAETRDFIDAYLKEMSKHTGNPTSSFHEENL ICSTLDLFFAGTETTSTTLRWALLYMALYPEIQEKVQAEIDRVIGQGQQPSTAA RESMPYTN AVIHEVQRMGNIIPLNVPREVTVDTTLAGYHLPKGT MILTNLTALH RDPTEWATPDTFNP DHFLENGQFKKREAFMPFSIGKRACLGEQLARTELFIF F TSLMQKFTFRPPNNEKLSLKFRMGITISPVSHRLCAVPQV
Research Area	Cardiovascular
Source	E.coli
Target Names	CYP2J2
Expression Region	33-502aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	61.3 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.