




Recombinant Human Probable ATP-dependent RNA helicase DDX5 (DDX5)

Product Code	CSB-EP006630HU
Relevance	Involved in the alternative regulation of pre-mRNA splicing; its RNA helicase activity is necessary for increasing tau exon 10 inclusion and occurs in a RBM4-dependent manner. Binds to the tau pre-mRNA in the stem-loop region downstream of exon 10. The rate of ATP hydrolysis is highly stimulated by single-stranded RNA. Involved in transcriptional regulation; the function is independent of the RNA helicase activity. Transcriptional coactivator for androgen receptor AR but probably not ESR1. Synergizes with DDX17 and SRA1 RNA to activate MYOD1 transcriptional activity and involved in skeletal muscle differentiation. Transcriptional coactivator for p53/TP53 and involved in p53/TP53 transcriptional response to DNA damage and p53/TP53-dependent apoptosis. Transcriptional coactivator for RUNX2 and involved in regulation of osteoblast differentiation. Acts as transcriptional repressor in a promoter-specific manner; the function probably involves association with histone deacetylases, such as HDAC1. As component of a large PER complex is involved in the inhibition of 3' transcriptional termination of circadian target genes such as PER1 and NR1D1 and the control of the circadian rhythms.
Abbreviation	Recombinant Human DDX5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17844
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSGYSSDRDRGRDRGFGAPRFGGSRAGPLSGKKFGNPGEKLVKKKWNLDEL LPKFEKNFYQEHPDLARRTAQEVETVRRSKEITVRGHNCPPVLFNFYEANFPA NVMDVIARQNFTEPTAIQAQGWVVALSGLDMVGVAQTGSGKTL SYLLPAIVHI NHQPFLERGDGPICLVLAPTRELAQQVQQVAAEYCRACRLKSTCIYGGAPKGP QIRDLERGVEICIA TPGR LIDFLECGKTNLRRTTYLV LDEADRMLDMGFEPQIRK IVDQIRPDRQTL MWSATWPKEVRQLAEDFLKDYIHINIGALELSANHNILQIVDV CHDVEKDEKLIRLMEEIMSEKENKTIVFVETKRRCEDELTRKMRRDGWPAMGIH GDKSQQERDWLNEFKH GKAPILIATDVASRGLDVEDVKFVINYDYPNSSSEDI HRIGRTARSTKTGTAYTFFTPNNIKQVSD LISVLREANQAINPKLLQLVEDRGSG RSRGRGGMKDDRRDRYSAGKRGGFNTFRDRENYDRGYSSLLKRDFGAKTQ NGVYSAANYTNGSFGSNFVSAGIQTSFRTGNPTGT YQNGYDSTQQYGSNVP NMHNGMNQQAYAYPATAAAPMIGYPMPTGYSQ
Research Area	Epigenetics and Nuclear Signaling




Source	E.coli
Target Names	DDX5
Protein Names	DEAD box protein 5 (RNA helicase p68 G17P1) (HELR) (HLR1)
Expression Region	1-614aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	74.1 kDa
Protein Length	Full Length


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP006630HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) DDX5.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP006630HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) DDX5.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.