




Recombinant Human GTP-binding protein Di-Ras3 (DIRAS3)

Product Code	CSB-EP006907HU
Abbreviation	Recombinant Human DIRAS3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95661
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MGNASFGSKEQKLLKRLRLLPALLILRAFKPHRKIRDYRVVVVGTTAGVGKSTLL HKWASGNFRHEYLP TIENTYCQLLGCSHGVL SLHITDSKSGDGNRALQRHVIA RGHAFVLVYSVTKKETLEELKAFYELICKIKGNNLHKFP IVLVGNKSDDTHREVA LNDGATCAMEWNC AFMEISAKTDVNVQELFHMLLNYKKKPTTGLQEPEKKSQ MPNTTEKLLDKCIIM
Research Area	others
Source	E.coli
Target Names	DIRAS3
Protein Names	Recommended name: GTP-binding protein Di-Ras3 Alternative name(s): Distinct subgroup of the Ras family member 3 Rho-related GTP-binding protein Rho1
Expression Region	1-229aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	52.8 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.