




Recombinant Mouse Desmoglein-3 (Dsg3), partial

Product Code	CSB-EP007205MOe1
Abbreviation	Recombinant Mouse Dsg3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O35902
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	EWVKFAKPCREREDNSRRNPIAKITSDFQKNQKITYRISGVGIDQPPFGIFVVD PNNGDINITAIVDREETPSFLITCRALNALGQDVERPLILTVKILDVNDNPPIFSQT IFKGEIEENSASNSLVMILNATDADEPNHMHNSKIAFKIVSQEPAGMSMFLISRNT GEVRTLTSSLDREQISSYHLVVGADNDGTGLSTQCECSIKIKDVNDNFPVLRE SQYSARIEENTLNAELLRFQVTDWDEEYTDNWLAVYFFTSGNEGNWFEIETDP RTNEGILKVVKALDYEQVQSMQFSIAVRNKAEFHQSVISQYRVQSTPVTIQVID VREGISFRPPSKTFTVQRGVSTNKLVG YILGT YQATDEDTGKAASSVRYVLGR NDGGLLVIDSKTAQIKFVKNIDRDSTFIVNKTISAEVLIDENTGKTSTGTIYVEVP SFNENCPSVVLEKKDICTSSPSVTL SVRTLDRGKYTG PYTVSLEEQLKLPVM WTITTLNATSALLQAQQQVSPGVYNVPVIVKDNQDGLCDTPESLTLTVQCDD RSMCRAPIPSREPNTYGESSWRLGP
Research Area	Others
Source	E.coli
Target Names	Dsg3
Expression Region	50-617aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	63.2 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.