




Recombinant Human Eukaryotic translation initiation factor 3 subunit E (EIF3E)

Product Code	CSB-EP007534HU
Relevance	omponent of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis . The eIF-3 complex associates with the 40S ribosome and facilitates the recruitment of eIF-1, eIF-1A, eIF-2:GTP:methionyl-tRNAi and eIF-5 to form the 43S pre-initiation complex (43S PIC). The eIF-3 complex stimulates mRNA recruitment to the 43S PIC and scanning of the mRNA for AUG recognition. The eIF-3 complex is also required for disassembly and recycling of post-termination ribosomal complexes and subsequently prevents premature joining of the 40S and 60S ribosomal subunits prior to initiation . The eIF-3 complex specifically targets and initiates translation of a subset of mRNAs involved in cell proliferation, including cell cycling, differentiation and apoptosis, and uses different modes of RNA stem-loop binding to exert either translational activation or repression . Required for nonsense-mediated mRNA decay (NMD); may act in conjunction with UPF2 to divert mRNAs from translation to the NMD pathway . May interact with MCM7 and EPAS1 and regulate the proteasome-mediated degradation of these proteins.
Abbreviation	Recombinant Human EIF3E protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P60228
Alias	Eukaryotic translation initiation factor 3 subunit 6 Viral integration site protein INT-6 homolog eIF-3 p48
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAEYDLTTRIAHFLDRHLVFPILLEFLSVKEIYNEKELLQGKLDLLSDTNMVDFA MDVYKNLYSDDIPHALREKRTTVVAQLKQLQAETEPIVKMFEDPETTRQMQST RDGRMLFDYLADKHGFRQEYLDTRYAKFQYECGNYSGAAYLYFFRVLP ATDRNALSSLWGKLASEILMQNWDAAEMDLRLKETIDNNSVSSPLQSLQQRT WLIHWSLFFVFNHPKGRDNIIDFLYQPQYLNAIQTMCPHILRYLTTAVITNKDV RKRRQVLKDLVKVIQQESYTYKDPITEFVECLYVNFDFDGAQKKLRECESVLV NDFFLVACLEDFIENARLFIFETFCRIHQCSINMLADKLNMTPEEAERWIVNLIR NARLDAKIDSKLGHVVMGNNVSPYQQVIEKTKSLSFRSQMLAMNIEKKLNQN SRSEAPNWATQDSGFY
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli



Target Names	EIF3E
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit E Short name= eIF3eAlternative name(s): Eukaryotic translation initiation factor 3 subunit 6 Viral integration site protein INT-6 homolog eIF-3 p48
Expression Region	1-445aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	79.1kDa
Protein Length	Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.