




Recombinant Rat Erythropoietin receptor (Epor), partial

Product Code	CSB-EP007744RA
Abbreviation	Recombinant Rat Epor protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q07303
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ASSPSLPDPKFESKAALLASRGSEELLCFTQRLDLVCFWEEAANS GMGFNY SFSYQLEGESRKSCRLHQAPTVRGSMRFWCSLPTADTSS FVPLELQVTEASG SPRYHRIIHINEVLLDAPAGLLARRAEEGSHVVLRLWLP PP GAPMTTHIRYEVD VSAGNRAGGTQRVEVLEGRTECVLSNLRGGTRYTF AVRARMAEPSFS GFWS AWSEPASLLTASDLDP
Research Area	Cardiovascular
Source	E.coli
Target Names	Epor
Expression Region	25-249aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	25.7 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.