




# Recombinant Human Protein FAM3B (FAM3B)

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP008227HU                                                                                                                                                                                                                                                                  |
| <b>Abbreviation</b>      | Recombinant Human FAM3B protein                                                                                                                                                                                                                                                 |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | P58499                                                                                                                                                                                                                                                                          |
| <b>Form</b>              | Liquid or Lyophilized powder                                                                                                                                                                                                                                                    |
| <b>Storage Buffer</b>    | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.                                                      |
| <b>Product Type</b>      | Recombinant Human Protein FAM3B(FAM3B)                                                                                                                                                                                                                                          |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                     |
| <b>Sequence</b>          | ELIPDAPLSSAAYSIRSIGERPVLKAPVPKRQKCDHWTPCPSDTYAYRLLSGGG<br>RSKYAKICFEDNLLMGEQLGNVARGINIAIVNYVTGNVTATRCFDMYEGDNSG<br>PMTKFIQSAAPKSLLFMVTYDDGSTRLNNDKNAIEALGSKEIRNMKFRSSWV<br>FIAAKGLELPSEIQREKINHSDAKNNRYSGWPAEIQIEGCIPKERS                                                       |
| <b>Research Area</b>     | Immunology                                                                                                                                                                                                                                                                      |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                          |
| <b>Target Names</b>      | FAM3B                                                                                                                                                                                                                                                                           |
| <b>Expression Region</b> | 30-235aa                                                                                                                                                                                                                                                                        |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal 10xHis-tagged and C-terminal Myc-tagged                                                                                                                                                                                                                              |
| <b>Mol. Weight</b>       | 30.3 kDa                                                                                                                                                                                                                                                                        |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                   |

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.