




Recombinant Human Low affinity immunoglobulin epsilon Fc receptor (FCER2), partial

Product Code	CSB-EP008534HU
Relevance	Low-affinity receptor for immunoglobulin E (IgE) and CR2/CD21. Has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).
Abbreviation	Recombinant Human FCER2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06734
Alias	BLAST-2;C-type lectin domain family 4 member JFc-epsilon-RIIImmunoglobulin E-binding factorLymphocyte IgE receptor; CD23
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	DTTQSLKQLEERAARNVSQVSKNLESHHGDQMAQKSQSTQISQELEELRAEQ QRLKSQDLELSWNLNGLQADLSSFKSQELNERNEASDLLERLREEVTKLRMEL QVSSGFVCNTCPEKWINFQRKCYFYGKGTKQVWHARYACDDMEGQLVSIHS PEEQDFLTKHASHTGSWIGLRNLDLKGFEIWDGSHVDYSNWAPGEPTSRSQ GEDCVMMRGSGRWNDAFCDRKLGAWVCDRLATCTPPASEGSAESMGPD SRPDPDGRLLPTPSAPLHS
Research Area	Immunology
Source	E.coli
Target Names	FCER2
Expression Region	48-321aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	47.0kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The gene encoding the Human FCER2 protein (48-321aa) creates recombinant plasmid, which is then transformed into e.coli cells. e.coli cells that can endure a specific antibiotic are selected and cultured under conditions that encourage the expression of the gene of interest. The protein features a N-terminal 6xHis-SUMO tag fusion. Following expression, affinity purification is employed to isolate and purify the recombinant Human FCER2 protein from the cell lysate. Denaturing SDS-PAGE is applied to resolve the resulting recombinant Human FCER2 protein, indicating a purity level exceeding 90%.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.