




Recombinant Human Farnesyl pyrophosphate synthase (FDPS)

Product Code	CSB-EP008563HU
Relevance	Key enzyme in isoprenoid biosynthesis which catalyzes the formation of farnesyl diphosphate (FPP), a precursor for several classes of essential metabolites including sterols, dolichols, carotenoids, and ubiquinones. FPP also serves as substrate for protein farnesylation and geranylgeranylation. Catalyzes the sequential condensation of isopentenyl pyrophosphate with the allylic pyrophosphates, dimethylallyl pyrophosphate, and then with the resultant geranylpyrophosphate to the ultimate product farnesyl pyrophosphate.
Abbreviation	Recombinant Human FDPS protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14324
Alias	(2E,6E)-farnesyl diphosphate synthase Dimethylallyltranstransferase (EC:2.5.1.1) Farnesyl diphosphate synthase Geranyltranstransferase
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MPLSRWLRSVGVFLLPAPYWAPRERWLGSLRRPSLVHGYPVLAWHSARCWC QAWTEEPRALCSSLRMNGDQNSDVYAQEKQDFVQHFSQIVRVLTEDMGHP EIGDAIARLKEVLEYNAIGGKYNRGLTVVVAFRELVEPRKQDADSLQRAWTVG WCVELLQAFFLVADDIMDSSLTRRGQICWYQKPGVGLDAINDANLLEACIYRLL KLYCREQPYLNLIELFLQSSYQTEIGQTLDLLTAPQGNVDLVRFTKRYKSIVK YKTAFYSFYLPAAAMYMAGIDGEKEHANAKKILLEMGEFFQIQDDYDLFGDP SVTGKIGTDIQDNKCSWLVVQCLQRATPEQYQILKENYQGQKEAEKVARVKALY EELDPAVFLQYEEDSYSHIMALIEQYAAPLPPAVFLGLARKIYKRRK
Research Area	Metabolism
Source	E.coli
Target Names	FDPS
Expression Region	1-419aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	64.3kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.