




Recombinant Mouse Fibroblast growth factor 20 (Fgf20)

Product Code	CSB-EP008626MOa0
Abbreviation	Recombinant Mouse Fgf20 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ESL9
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MAPLTEVGAFLLGGLEGLSQVVGSHFLLPPAGERPPLLGERRGALERGARGGP GSVELAHLHGILRRRQLYCRTGFHLQILPDGTVQGTRQDHSFLFGILEFISVAVG LVSIRGVDSGLYLG MNDKGELYGSEKLTSECIFREQFEENWYNTYSSNIYKHG DTGRRYFVALNKDGT PRDGARSKRHQKFTHFLPRPVDPERVPELYKDLLMYT
Research Area	Cancer
Source	E.coli
Target Names	Fgf20
Expression Region	1-211aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	27.7 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.