



Recombinant Human Fibroblast growth factor 3 (FGF3)

Product Code	CSB-EP008630HU
Abbreviation	Recombinant Human FGF3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11487
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Fibroblast growth factor 3(FGF3)
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	AAGPGARLRRDAGGRGGVYEHLLGGAPRRRKLYCATKYHLQLHPSGRVNGSL ENSAYSILEITAVEVGIVAIRGLFSGRYLAMNKRGRLYASEHYSAECECFVERIHE LGYN TYASRLYRTVSSTPGARRQPSAERLWYVSVNGKGRPRRGFKTRRTQK SSLFLPRVLDHRDHEMVRQLQSGLP RPPGKGVQPRRRRQKQSPDNLEPSHV QASRLGSQLEASAH
Research Area	Others
Source	E.coli
Target Names	FGF3
Expression Region	18-239aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	32.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The recombinant Human FGF3 protein is produced by manipulating the expression of specific genes in e.coli cells. A DNA sequence encoding the Human FGF3 protein (18-239aa) is fused to an expression vector and subsequently transformed into e.coli cells. The cells that successfully integrate the vector are selected and cultured to induce the expression of the desired protein. A N-terminal 10xHis tag and C-terminal Myc tag is attached to the protein. The recombinant Human FGF3 protein is then purified using affinity purification, achieving a purity level of over 85%, as confirmed by SDS-PAGE analysis.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.