




Recombinant Human Fragile X mental retardation syndrome-related protein 1 (FXR1)

Product Code	CSB-EP009087HU(A4)e5
Abbreviation	Recombinant Human FXR1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51114
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Fragile X mental retardation syndrome-related protein 1 (FXR1)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAELTVEVRGSNGAFYKGFIKDVHEDSLTVVFENNWQPERQVPFNEVRLPPP PDIKKEISEGDEVEVYSRANDQEPGWWLAKVRMMKGEFYVIEYAACDATYN EIVTFERLRPVNQNKTVKKNTFFKCTVDVPEDLREACANENAHKDFKKAVGAC RIFYHPETTQLMILSASEATVKRVNLSDMHLRSIRTKLMLMSRNEEATKHLECT KQLAAAFHEEFVVREDLMGLAIGTHGSNIQQARKVPGVTAIELDEDTGTFRIYG ESADAVKKARGFLEFVEDFIQVPRNLVGKVIQKNGKVIQEIQVDSGVVRVRIEG DNENKLPREDGMVPFVFGTKESIGNVQVLLEYHIAYLKEVEQLRMERLQIDE QLRQIGSRYSRGRGRGRRPNYTSYGVTNSELSNPSETESERKDELSDSL AGEDDRDSRHRQDSRRRPGGRGRSVSGGRGRGGPRGGKSSISSVLKDPDS NPYSLLDNTESDQTADTDASESHHSTNRRRRSRRRRTDEDAVLMDGMTESD TASVNENGLTVADYISRAESQSRQRNLPRETLAKNKKEMAKDVIEEHGPSEK AINGPTSASGDDISKLRTPGEEKINTLKEENTQEAAVLNGVS
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	FXR1
Expression Region	1-621aa
Tag Info	N-terminal GFP-tagged and C-terminal 6xHis-tagged
Mol. Weight	97.7 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.