




Recombinant Human Phospholipid hydroperoxide glutathione peroxidase GPX4 (GPX4) (U73S), partial

Product Code	CSB-EP009869HU
Abbreviation	Recombinant Human GPX4 protein (U73S), partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36969
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	CASRDDWRCARSMHEFSAKDIDGHMVNLDKYRGFVCIVTNVASQSGKTEVNY TQLVDLHARYAECGLRILAFPCNQFGKQEPGSNEEIKEFAAGYNVKFDMFSKI CVNGDDAHPLWKWMKIQPKGKGILGNAIKWNFTKFLIDKNGCVVKRYGPMEE PLVIEKDLPHYF
Research Area	Cancer
Source	E.coli
Target Names	GPX4
Expression Region	29-197aa(U73S)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	23.4 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The recombinant human GPX4 protein with an N-terminal 6xHis tag is expressed using an E. coli system. The GPX4 gene fragment [29-197aa (U73S)] is co-cloned into an expression plasmid with the tag gene. The recombinant plasmid is introduced into E. coli cells, which are cultured for protein expression. Protein expression is induced by IPTG, and the cells are lysed to release the recombinant GPX4 protein. Protein purification is achieved using Ni-NTA affinity chromatography, where the His tag binds to nickel ions, enabling the selective purification of GPX4. The eluted GPX4 protein is analyzed by SDS-PAGE, revealing a purity of greater than 90%.

Human GPX4 is a crucial selenoenzyme that plays a significant role in cellular antioxidant defense mechanisms, particularly in the reduction of lipid hydroperoxides and the maintenance of redox homeostasis. Its unique function as the only known enzyme capable of reducing phospholipid hydroperoxides distinguishes it from other glutathione peroxidases, which primarily target hydrogen peroxide and organic hydroperoxides [1][2]. The enzyme's activity is critically dependent on a selenocysteine residue, which is essential for its catalytic function [3].

GPX4 is involved in regulating ferroptosis [1][4]. The enzyme's ability to mitigate oxidative stress is vital for cell survival, particularly in tissues with high oxidative demands, such as the brain and the gastrointestinal tract [5]. GPX4's role extends beyond mere antioxidant activity. It also participates in signaling pathways that affect cell proliferation and apoptosis. In cancer biology, GPX4 has been identified as a potential therapeutic target, as its inhibition can enhance the efficacy of chemotherapeutic agents like cisplatin by promoting ferroptosis cell death in tumor cells [6]. This dual role of GPX4 in both protecting against oxidative damage and facilitating cell death in cancer cells highlights its importance in both health and disease.

References:

- [1] W. Yang, R. Sriramaratnam, M. Welsch, K. Shimada, R. Skouta, V. Viswanathan, et al. Regulation of ferroptotic cancer cell death by gpx4, *Cell*, vol. 156, no. 1-2, p. 317-331, 2014. <https://doi.org/10.1016/j.cell.2013.12.010>
- [2] J. Angeli, M. Schneider, B. Proneth, Y. Tyurina, V. Tyurin, V. Hammond, et al. Inactivation of the ferroptosis regulator gpx4 triggers acute renal failure in mice, *Nature Cell Biology*, vol. 16, no. 12, p. 1180-1191, 2014. <https://doi.org/10.1038/ncb3064>



[3] S. Zhu, Q. Zhang, X. Sun, H. Zeh, M. Lotze, R. Kang, et al. Hspa5 regulates ferroptotic cell death in cancer cells, *Cancer Research*, vol. 77, no. 8, p. 2064-2077, 2017. <https://doi.org/10.1158/0008-5472.can-16-1979>

[4] M. Gaschler, A. Andia, H. Liu, J. Csuka, B. Hurlocker, C. Vaiana et al., Fno2 initiates ferroptosis through gpx4 inactivation and iron oxidation, *Nature Chemical Biology*, vol. 14, no. 5, p. 507-515, 2018. <https://doi.org/10.1038/s41589-018-0031-6>

[5] B. Speckmann, H. Bidmon, A. Pinto, M. Anlauf, H. Sies, & H. Steinbrenner, Induction of glutathione peroxidase 4 expression during enterocytic cell differentiation, *Journal of Biological Chemistry*, vol. 286, no. 12, p. 10764-10772, 2011. <https://doi.org/10.1074/jbc.m110.216028>

[6] X. Zhang, S. Sui, L. Wang, H. Li, L. Zhang, S. Xu, et al. Inhibition of tumor propellant glutathione peroxidase 4 induces ferroptosis in cancer cells and enhances anticancer effect of cisplatin, *Journal of Cellular Physiology*, vol. 235, no. 4, p. 3425-3437, 2019. <https://doi.org/10.1002/jcp.29232>

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.