




Recombinant Human Histone H3.1 (HIST1H3A), partial

Product Code	CSB-EP010418HU1
Abbreviation	Recombinant Human HIST1H3A protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P68431
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ARTKQTARKSTGGKAPRKQLATKAARKSAPATGGVKKPHRYRPGTVALREIR RYQKSTELLIRKLFPQRLVREIAQDFKTDLRFQSSAVMALQEACEAYLVGLFED TNLCAIHAKRVTIMPKDIQLARRIRGER
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	HIST1H3A
Expression Region	2-135aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	19.3 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.