




Recombinant Human Heparan-sulfate 6-O-sulfotransferase 1 (HS6ST1), partial

Product Code	CSB-EP010756HU1
Abbreviation	Recombinant Human HS6ST1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60243
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Heparan-sulfate 6-O-sulfotransferase 1 (HS6ST1), partial
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	PGLSLGAPGGRAPDDLDLFTPDPHYEKYYFPVRELERSLRFDMMKGDDVIV FLHIQKTGGTTFGRLHVQNVRLVPCDCRPGQKKCTCYRPNRRETWLFSSRS TGWSCGLHADWTELTNCVPGVLDRRDSAAALRTPRKFFYYITLLRDPVSRYLSE WRHVQRGATWKTSLHMC DGRTPTEELPPCYEGTDWSGCTLQEFMDCPYN LANNRQVRMLADLSLVGCYNLSFIPEGKRAQLLLES AKKNLRGMAFFGLTEFQ RKTQYLFERTFNLKFIRPFMQYNSTRAGGVEVDEDTIRRIEELNDLDMQLYDYA KDLFQQRYQYKRQLERREQRLRSREERLLHRAKEALPREDADEPGRVPTEDY MSHIEKW
Research Area	Neuroscience
Source	E.coli
Target Names	HS6ST1
Expression Region	38-411aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	51.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.