




Recombinant Mouse Ketohexokinase (Khk)

Product Code	CSB-EP012157MO
Abbreviation	Recombinant Mouse Khk protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P97328
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MEEKQILCVGLVVLDIINVVDKYPEEDTDRRCLSQRWQRGGNASNSCTVLSLL GARCAFMGSLAPGHVADFLVADFRQRGVDVSQVTWQSQGDTPCSCCIVNNS NGSRTIILYDTNLPDVSAKDFEKVDLTRFKWIHIEGRNASEQVKMLQRIEEHNA KQPLPQKVRVSVEIEKPREELFQLFSYGEVVFVSKDVAKHLGFQPAVEALRGL YSRVKKGATLVCAWAEEGADALGPDGQLLHSDAFPPPRVVDTLGAGDTFNAS VIFSLSKGNSMQEALRFGCQVAGKKCGLQGFDGIV
Research Area	Cancer
Source	E.coli
Target Names	Khk
Protein Names	Recommended name: Ketohexokinase EC= 2.7.1.3 Alternative name(s): Hepatic fructokinase
Expression Region	1-298aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	36.4kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.