




Recombinant Human Keratin,type II cytoskeletal 8 (KRT8)

Product Code	CSB-EP012574HU(A4)
Abbreviation	Recombinant Human KRT8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05787
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSIRVTQKSYKVSTSGPRAFSSRSYTS GP GSRISSSSF SRVGSSNFRGGLGG GYGGASGMGGITAVTVNQSLLSPLVLEVD PNIQAVRTQEKEQIKTLNNKFASFI DKVRFLEQQNKMLETKWSLLQQQKTARSNMDNMFESYINNLRRLQLET LGQEK LKLEAELGNMQGLVEDFKNKYEDEINKRTEMENEFVLIKKDVDEAYMNKVELE SRLEGLTDEINFLRQLYEEEIRELQSQISDTSVVL SMDNSRSLDMSIIAEVKAQ YEDIANRSRAEAESMYQIKYEELQSLAGKHGDDLRR TKTEISEMNRNISRLQAE IEGLKGQRASLEAAIADAEQRGELAIKDANAKLSELEAALQRAKQDMARQLRE YQELMNVKLALDIEIATYRKLLGEESRLES GMQNMSIHTKTTSGYAGGLSSAY GGLTSPGLSYSLGSSFGSGAGSSSF SRTSSSRVVVKIETRDGKLVSESSDV LPK
Research Area	Signal Transduction
Source	E.coli
Target Names	KRT8
Expression Region	1-483aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	57.8 kDa
Protein Length	Full length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.