




Recombinant Mouse Lipoprotein lipase (Lpl)

Product Code	CSB-EP013065MOe1
Abbreviation	Recombinant Mouse Lpl protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11152
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ADAGRDFSDIESKFALRTPEDTAEDTCHLIPGLADSVSNCHFHNHSSKTFVVIHG WTVTGMYESWVPKLVAALYKREPDSNVIVVDWLYRAQQHYVVSAGYTKLVG NDVARFINWMEEEFNYPLDNVHLLGYSLGAHAAGVAGSLTNKKVNRITGLDPA GPNFEYAEAPSRLSPDDADFVDVLHTFTRGSPGRSIGIQKPVGHVDIYPNGGT FQPGCNIGEAIRVIAERGLGDVDQLVKCSHERSIHLFIDSLLNEENPSKAYRCN SKEAFEKGLCLSCRKNRCNNLGYEINKVRAKRSSKMYLKTRSQMPYKVFHYQ VKIHFSGTEDGKQHNQAFEISLYGTVAESENIPFTLPEVSTNKTYSFLIYTEVDIG ELLMMKLKWISDSYFSWPDWWSSPSFVIERIRVKAGETQKKVIFCAREKVS HLQKGKDSAVFVKCHDKSLKKSG
Research Area	Cancer
Source	E.coli
Target Names	Lpl
Expression Region	28-474aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	50.3 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.