




Recombinant Mouse Dual specificity mitogen-activated protein kinase kinase 4 (Map2k4)

Product Code	CSB-EP013413MO
Relevance	Dual specificity protein kinase which acts as an essential component of the MAP kinase signal transduction pathway. Essential component of the stress-activated protein kinase/c-Jun N-terminal kinase (SAP/JNK) signaling pathway. With MAP2K7/MKK7, is the one of the only known kinase to directly activate the stress-activated protein kinase/c-Jun N-terminal kinases MAPK8/JNK1, MAPK9/JNK2 and MAPK10/JNK3. MAP2K4/MKK4 and MAP2K7/MKK7 both activate the JNKs by phosphorylation, but they differ in their preference for the phosphorylation site in the Thr-Pro-Tyr motif. MAP2K4 shows preference for phosphorylation of the Tyr residue and MAP2K7/MKK7 for the Thr residue. The phosphorylation of the Thr residue by MAP2K7/MKK7 ses to be the prerequisite for JNK activation at least in response to proinflammatory cytokines, while other stimuli activate both MAP2K4/MKK4 and MAP2K7/MKK7 which synergistically phosphorylate JNKs. MAP2K4 is required for maintaining peripheral lymphoid homeostasis. The MKK/JNK signaling pathway is also involved in mitochondrial death signaling pathway, including the release cytochrome c, leading to apoptosis. Whereas MAP2K7/MKK7 exclusively activates JNKs, MAP2K4/MKK4 additionally activates the p38 MAPKs MAPK11, MAPK12, MAPK13 and MAPK14.
Abbreviation	Recombinant Mouse Map2k4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47809
Alias	C-JUN N-terminal kinase kinase 1 ;JNK kinase 1 ;JNKK 1JNK-activating kinase 1MAPK/ERK kinase 4 ;MEK 4SAPK/ERK kinase 1 ;SEK1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	AAPSPSGGGGSGGGGGTPGPIGPPASGHPAVSSMQGKRKALKLNFANPPVK STARFTLNPNTTGTVQNPHIERLRTHSIESSGKLKISPEQHWDFTAEDLKDLGEI GRGAYGSVNKMVHKPSGQIMAVKRIRSTVDEKEQKQLLMDLDVVMRSSDCPY IVQFYGALFREGDCWICMELMSTSFDFKYVYVSVLDDVIPEEILGKITLATVKA LNHLKENLKIHRDIKPSNILLDRSGNIKLCDFGISGQLVDSIAKTRDAGCRPYMA PERIDPSASRQGYDVRSVDVWSLGITLYELATGRFPYPKWNSVFDQLTQVVKG DPPQLSNSEEREFSPSFINFVNLCLTKDESKRPKYKELLKHPFILMYEERTVEV ACYVCKILDQMPATPSSPMYVD
Research Area	Others



Source	E.coli
Target Names	Map2k4
Expression Region	2-397a
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	60.0kDa
Protein Length	Full Length of Mature Protein

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.