




Recombinant Rat Beta-nerve growth factor (Ngf), partial

Product Code	CSB-EP015779RA1
Relevance	Nerve growth factor is important for the development and maintenance of the sympathetic and sensory nervous systems. Extracellular domain ligand for the NTRK1 and NGFR receptors, activates cellular signaling cascades through those receptor tyrosine kinase to regulate neuronal proliferation, differentiation and survival. Inhibits metalloproteinase dependent proteolysis of platelet glycoprotein VI.
Abbreviation	Recombinant Rat Ngf protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25427
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Alias	Beta-NGF
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	THPVFHMGEFSVCDSVSVWVGDKTTATDIKGKEVTVLGEVNINNSVFKQYFFE TKCRAPNPVESGCRGIDSKHWNSYCTTTHTFVKALTDDKQAAWRFIRIDTAC VCVLSRKAARRG
Research Area	Others
Source	E.coli
Target Names	Ngf
Expression Region	124-241aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.1 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.