




Recombinant Human Nucleolar protein 3 (NOL3)

Product Code	CSB-EP015921HU
Relevance	Isoform 1: May be involved in RNA splicing.
Abbreviation	Recombinant Human NOL3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60936
Alias	Apoptosis repressor with CARD1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MGNAQERPSETIDRERKRLVETLQADSGLLLDALLARGVLTGPEYEALDALPD AERRVRLLLLLVQKGGEAACQELLRCAQRTAGAPDPAWDWQHVGPGYRDR SYDPPCPGHWTPPEAPGSGTTCPLPRASDPDEAGGPEGSEAVQSGTPPEEPE PELEAEASKEAEPEPEPEPELEPEAEAEPEPEPEPEPDPEPEPDFEERDESED S
Research Area	Apoptosis
Source	E.coli
Target Names	NOL3
Expression Region	1-208aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	38.6kDa
Protein Length	Full Length of Isoform 2

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.