




Recombinant Human Puromycin-sensitive aminopeptidase (NPEPPS), partial

Product Code	CSB-EP015980HU
Abbreviation	Recombinant Human NPEPPS protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55786
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	RLPADVSPINYSLCLKPDLLDFTFEGKLEAAAQVRQATNQIVMNCADIDIITASY APEGDEEIHATGFNYQNEDEKVTLSFPSTLQTGTGTLKIDFVGELNDKMKGFY RSKYTTSPSGEVRYAAVTQFEATDARRAFPCWDEPAIKATFDISLVVPKDRVAL SNMNVIDRKPYPDDENLVEVKFARTPVMSTYLVAFFVVGGEYDFVETRSDKGVC VRVYTPVGKAEQGGKFALEVAAKTLPFYKDYFNVPYPLPKIDLIAIDFAAGAME NWGLVTYRETALLIDPKNSCSSSRQWVALVVGHELAHQWFGNLVTMEWWTH LWLNEGFASWIEYLCVDHCFPEYDIWTQFVSADYTRAQELDALDNSHPIEVSV GHPSEVDEIFDAISYSGASVIRMLHDYIGDKDFKKGMNMYLTKFQQKNAATE DLWESLENASGKPIAAVMNTWTKQMGFP
Research Area	Neuroscience
Source	E.coli
Target Names	NPEPPS
Protein Names	Recommended name: Puromycin-sensitive aminopeptidase Short name= PSA EC= 3.4.11.14Alternative name(s): Cytosol alanyl aminopeptidase Short name= AAP-S
Expression Region	53-504aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	57.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.