




Recombinant Human Nuclear receptor subfamily 2 group F member 6 (NR2F6)

Product Code	CSB-EP016058HU
Relevance	Transcription factor predominantly involved in transcriptional repression. Binds to promoter/enhancer response elements that contain the imperfect 5'-AGGTCA-3' direct or inverted repeats with various spacings which are also recognized by other nuclear hormone receptors. Involved in modulation of hormonal responses. Represses transcriptional activity of the lutropin-choriogonadotropic hormone receptor/LHCGR gene, the renin/REN gene and the oxytocin-neurophysin/OXT gene. Represses the triiodothyronine-dependent and -independent transcriptional activity of the thyroid hormone receptor gene in a cell type-specific manner. The corepressing function towards thyroid hormone receptor beta/THRB involves at least in part the inhibition of THRB binding to triiodothyronine response elements (TREs) by NR2F6. Inhibits NFATC transcription factor DNA binding and subsequently its transcriptional activity. Acts as transcriptional repressor of IL-17 expression in Th-17 differentiated CD4+ T cells and may be involved in induction and/or maintenance of peripheral immunological tolerance and autoimmunity. Involved in development of forebrain circadian clock; is required early in the development of the locus coeruleus (LC).
Abbreviation	Recombinant Human NR2F6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10588
Alias	V-erbA-related protein 2 Short name: EAR-2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAMVTGGWGGPGGDTNGVDKAGGYPRAAEDDSASPPGAASDAEPGDEERP GLQVDCVVC GDKSSGKH YGVFTCEGCKSFFKRSIRRNL SYTCRSNRDCQIDQ HHRNQCYCRLKKCFRVGMRKEAVQRGRIPHS LPGA VAASSGSPPGSALAA VASGGDLFP GQPVSELIAQLLRAEPYPAAAGRFGAGGGAAGAVLGIDNVCELA ARLLFSTVEWARHAPFFPELPVADQVALLRLSWSELFVLNAAQAALPLHTAPLL AAAGLHAAPMAAERAVAFMDQVRAFQE QVDKLGRLQVDSA EYGCLKAIALFT PDACGLSDPAHVESLQEKAQVALTEYVRAQYPSQPQRFGRLLLR L PALRAVP ASLISQLFFMRLVGKTP IETLIRDMLLSGSTFNWPY GSGQ
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli




Target Names	NR2F6
Expression Region	1-404aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	59.0kDa
Protein Length	Full Length


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP016058HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) NR2F6.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP016058HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) NR2F6.

Description

CUSABIO team inserts the gene coding for the Human NR2F6 protein (1-404aa) into a plasmid vector to form recombinant plasmid, which is then introduced into e.coli cells. e.coli cells demonstrating successful uptake of the recombinant plasmid are selected based on their ability to survive in the presence of a specific antibiotic. The positive e.coli cells are cultured under conditions that promote the expression of the gene of interest. A N-terminal 6xHis-SUMO tag is linked to the protein. Following expression, affinity purification is used to isolate and purify the recombinant Human NR2F6 protein from the cell lysate. Denaturing SDS-PAGE is then applied to resolve the resulting recombinant Human NR2F6 protein, demonstrating a purity greater than 90%.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.