




Recombinant Human Opioid growth factor receptor (OGFR)

Product Code	CSB-EP016311HU
Abbreviation	Recombinant Human OGFR protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZT2
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MDDPDCDSTWEEDEEDAEDAEDDCEDGEEAAGARDADAGDEDEESEEPRARPSSSFQSRMTGSRNWRATRDMDCRYRHNYPDLVERDCNGDTPNLSFYRNEIRFLPNGCFIEDILQNWTDNYDLLEDNHSYIQWLFPLREPGVNVHAKPLTLREV EVFKSSQEIQRRLVRAYELMLGFYGIRLEDRTGTGTVGRAQNYQKRFQNLNWR SHNNLRITRILKSLGELGLEHFQAPLVRFFLEETLVRRELPGVRQSALDYFMFA VRCRHQRRQLVHFAWEHFRPRCKFVWGPQDKLRRFKPSSLPHPLEGSRKVE EEGSPGDPDHEASTQGRTCGPEHSGGGGRVDEGPQPRSVEPQDAGPLERS QGDEAGGHGEDRPEPLSPKESKKRKLELSRREQPPTEPGPQSASEVEKIALN LEGALSQGSRLRTGTQEVGQDQPGAEVQPCRQPLGARVADKVRKRRKVDEG AGDSAAVASGGAQTLALAGSPAPSGHPKAGHSENGVEEDTEGRTGPKEGTP GSPSETPGPSPAGPAGDEPAESPSETPGPRPAGPAGDEPAESPSETPGPRP AGPAGDEPAESPSETPGPSPAGPTRDEPAESPSETPGPRPAGPAGDEPAESP SETPGPRPAGPAGDEPAESPSETPGPSPAGPTRDEPAKAGEAAELQDAEVES SAKSGKP
Research Area	Cell Biology
Source	E.coli
Target Names	OGFR
Protein Names	Recommended name: Opioid growth factor receptor Short name= OGFR Alternative name(s): Protein 7-60 Zeta-type opioid receptor
Expression Region	1-677aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	80.8 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.