



Recombinant Human DNA polymerase theta (POLQ), partial

Product Code	CSB-EP018321HU9a0
Abbreviation	Recombinant Human POLQ protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75417
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Sodium-dependent phosphate transport protein 2B (SLC34A2), partial
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SSSESLSIIDVASDQNLFQTFIKWRCKKRFSISLACEKIRSLTSSKTATIGSRFK QASSPQEIPRDDGFPKGCDDTLVVGAVCWGGGRDAYYFSLQKEQKHSEISA SLVPPSLDPSLTLKDRMWYLQSLRKESDKECSVVIYDFIQSYKILLSCGISLE QSYEDPKVACWLLDPDSQEPTLHSIVTSFLPHELPLLEGMETSQGIQSLGLNA GSEHSGRYRASVESILIFNSMNQLNSLLQKENLQDVFRKVEMPSQYCLALLEL NGIGFSTAECESQKHIMQAKLDAIETQAYQLAGHSFSFTSSDDIAEVLFLLELKL PNREMKNQGSKKTGSTRRGIDNGRKLRLGRQFSTSKDVLNKLKALHPLPLGLI LEWRRITNAITKVVFPLQREKCLNPFLGMERIYPVSQSHTATGRITFTEPNIQNV PRDFEIKMPTLVGESPPSQAVGKGLLPMGRGKYKKGFSVNPRCQAQMEERA ADRGMPFSISMRHAFVFPFGGSILAADYSQLELRILAHLSHDRRLIQVLNTGAD VFRSIAAEWKMIPEESVGDDLRRQQAQKICYGIYGMGAKSLGEQMGIKENDAA CYIDSFKSRYTGINQFMTETVKNCKRDGFVQTLGRRRYLPGIKDNNPYRKAHA ERQAINIVQGSAAIVKIATVNIQKQLETFTFKSHGHREGMLQSDQTGLSR KRKLQGMFCPIRGFFILQLHDELLYEVAEEDVVQVAQIVKNEMESAVKLSVKL KVKVKIGASWGELKDFDV
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	POLQ
Expression Region	1820-2590aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info N-terminal 6xHis-tagged

Mol. Weight 90.6 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.