




Recombinant Human Melanoma antigen preferentially expressed in tumors (PRAME)

Product Code	CSB-EP018603HU
Relevance	Functions as a transcriptional repressor, inhibiting the signaling of retinoic acid through the retinoic acid receptors RARA, RARB and RARG. Prevents retinoic acid-induced cell proliferation arrest, differentiation and apoptosis.
Abbreviation	Recombinant Human PRAME protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P78395
Alias	Opa-interacting protein 4 ;OIP-4Preferentially expressed antigen of melanoma
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MERRRLWGSIQSRYISMSVWTSPRRLVELAGQSLLKDEALAIAALELLPRELFP PLFMAAFDGRHSQTLKAMVQAWPFTCLPLGVLMKGQHLHLETFKAVLDGLDV LLAQEVRPRRWKLQVLDLRKNSHQDFWTVWSGNRASLYSFPEPEAAQPMTK KRKVDGLSTEAEQPFIPVEVLVDLFLKEGACDELFSYLIEKVKRKKNVLRCLCK KLKIFAMPMQDIKMILKMVQLDSIEDLEVTCTWKLPTLAKFSPYLGQMINLRRL LSHIHASSYISPEKEEQYIAQFTSQFLSLQCLQALYVDSLFFLRGRDLQLLRHVM NPLETLSITNCRLSEGDVMHLSQSPSVSQLSVLSLSGVMLTDVSPEPLQALLER ASATLQDLVFDECGITDDQLLALLPSLSHCSQLTTLSFYGNSISISALQSLLQHLL GLSNLTHVLYPVPLESYEDIHGTLHLERLAYLHARLRELLCELGRPSMVWLSAN PCPHCGDRTFYDPEPILCPCFMPN
Research Area	Apoptosis
Source	E.coli
Target Names	PRAME
Expression Region	1-509aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	61.9kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

Discover the power of our high-quality Recombinant Human MBL2 protein, engineered to provide optimal performance in immunology research. Expressed in *E. coli*, this protein represents the full length of the mature protein (2-115aa) and features an N-terminal 6xHis-tag for seamless purification. With a purity greater than 90% as determined by SDS-PAGE, our MBL2 protein ensures dependable and consistent results in your experimental endeavors.

Investigate the diverse roles of MBL2 in innate immunity, inflammation, and host-pathogen interactions using this exceptional recombinant protein. Available in both liquid and lyophilized powder forms, our Recombinant Human MBL2 protein is a valuable asset for enhancing your research capabilities.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.