




Recombinant Human Peroxiredoxin-1 (PRDX1)

Product Code	CSB-EP018653HU
Relevance	Involved in redox regulation of the cell. Reduces peroxides with reducing equivalents provided through the thioredoxin system but not from glutaredoxin. May play an important role in eliminating peroxides generated during metabolism. Might participate in the signaling cascades of growth factors and tumor necrosis factor-alpha by regulating the intracellular concentrations of H ₂ O ₂ . Reduces an intramolecular disulfide bond in GTPD5 that gates the ability to GTPD5 to drive postmitotic motor neuron differentiation
Abbreviation	Recombinant Human PRDX1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q06830
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSSGNAKIGH PAPNFKATAVMPDGQFKDISLSDYKGKYVVFFFYPLDFTFVCP TEIIAFSDRAEEFKLNCQVIGASVDSHFCHLAWVNTPKKQGGLGPMNIPLVSD PKRTIAQDYGV LKADEGISFRGLFIIDDKGILRQITVNDLPVGRSVDETLRLVQAF QFTDKHGEVCPAGWKPGSDTIKPDVQKSKEYFSKQK
Research Area	Cardiovascular
Source	E.coli
Target Names	PRDX1
Protein Names	Natural killer cell-enhancing factor A Proliferation-associated gene protein Thioredoxin peroxidase 2 Thioredoxin-dependent peroxide reductase 2 PAGA, PAGB, TDPX2
Expression Region	1-199aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	27.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.