




Recombinant Human 60S ribosomal protein L9 (RPL9)

Product Code	CSB-EP020326HU
Abbreviation	Recombinant Human RPL9 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32969
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MKTILSNQTVDIPENV DITLKGRTVIVKGP RGT LRRDFNHINVELSLLGKKKKRL RVDKWWGNRKELATVRTICSHVQNMIGVTLGFRYKMRSVYAHFPINVVIQEN GSLVEIRNFLGEKYIRRV RMRPGVACSVSQAQKDELILEGNDIELVSNSAALIQ QATTVKNKDIRKFLDGIYVSEKGT VQQADE
Research Area	others
Source	E.coli
Target Names	RPL9
Protein Names	Recommended name: 60S ribosomal protein L9
Expression Region	1-192aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	48.9kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.