




Recombinant Human 60S acidic ribosomal protein P0 (RPLP0)

Product Code	CSB-EP020336HU
Relevance	Ribosomal protein P0 is the functional equivalent of E.coli protein L10.
Abbreviation	Recombinant Human RPLP0 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05388
Alias	60S ribosomal protein L10E
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MPREDRATWKSNYFLKIIQLDDYPKCFIVGADNVGSKQMQQIRMSLRGKAVV LMGKNTMMRKAIRGHLENNPALEKLLPHIRGNVGFVFTKEDLTEIRDMLANKV PAAARAGAIAPCEVTVPAQNTGLGPEKTSFFQALGITTKISRGITIEILSDVQLIKT GDKVGASEATLLNMLNISPFSGFLVIQVQVFDNGSIYNPEVLDITEETLHSRFLEG VRNVASVCLQIGYPTVASVPHSIINGYKRVLALSVETDYTFPLAEKVKAFLADPS AFVAAAPVAAATTAAPAAAAAPAKVEAKEESESEDEDMGFGLFD
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	RPLP0
Expression Region	1-317aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	50.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The production of the recombinant Human RPLP0 protein begins with the creation of the recombinant plasmid, which is synthesized by inserting the gene encoding the Human RPLP0 protein (1-317aa) into a plasmid vector. The recombinant plasmid is introduced into e.coli cells. e.coli cells that can survive in the presence of a specific antibiotic are selected and then cultured under conditions conducive to the expression of the gene of interest. The protein is equipped with a N-terminal 6xHis-SUMO tag. Following expression, the recombinant RPLP0 protein is isolated and purified from the cell lysate using affinity purification. Denaturing SDS-PAGE is then employed to resolve the resulting recombinant Human RPLP0 protein, demonstrating a purity exceeding 90%.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.