




Recombinant Mouse Transcription factor SOX-1(Sox1)

Product Code	CSB-EP022417MO
Abbreviation	Recombinant Mouse Sox1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53783
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Transcription factor SOX-1(Sox1)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	<p> MYSMMMETDLHSPGGAQAPTNLSPAGAGGGGGGGGGGGGGGGGGTGANQD RVKRPMNAFMVWSRGQRRKMAQENPKMHNSEISKRLGAEWKVMSEAEKRP FIDEAKRLRALHMKEHPDYKYRPRRKTLLKDKYSLAGLLAAGAGGGGAA VAMGVGVGVGAAAVGQRLESPGGAAGGGYAHVNGWANGAYPGSVAAAAAA AAMMQEAQLAYGQHPGAGGAHPHAHPAHPHPHPAHPHPNPQPMHRYDM GALQYSPISNSQGYMSASPSGYGGIPYGAAGAGAGAGAGAGAGAGAGAGAG AAASSGALGALGSLVKSEPSGSPAPAHSRAPCPGDLREMISMYPAGEGGD PAAAAAAAAQSRHLHSLPQHYQGAGAGVNGTVPLTHI </p>
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Sox1
Expression Region	1-391aa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	46.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.