




Recombinant Human Signal recognition particle subunit SRP72 (SRP72)

Product Code	CSB-EP022677HU
Abbreviation	Recombinant Human SRP72 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O76094
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ASGGSGGVSVPALWSEVNRYGQNGDFTRALKTVNKLQINKDDVTALHCKVV CLIQNGSFKEALNVINTHTKVLANNLSFEKAYCEYRLNRIENALKTIESANQQT DKLKELYGQVLYRLERYDECLAVYRDLVRNSQDDYDEERKTNLSAVVAAQSN WEKVVPENLGLQEGTHELCYNTACALIGQGQLNQAMKILQKAEDLCRRSLSE DTDGTEEDPQAEIAIHGQMAYILQLQGRTEEALQLYNQIIKLKPTDVGLLAVIA NNIITINKDQNVFDSKKKVKLTNAEGVEFKLSKKQLQAIEFNKALLAMYTNAEQ CRKISASLQSQSPEHLLPVLIAAQLCREKQHTKAIELLQEFSDQHPENAAEIKL TMAQLKISQGNISKACLILRSIEELKHKPGMVSAVMTMYSHEEDIDSAIEVFTQAI QWYQNHQPKSPAHLSLIREAANFKLKYGRKKEAISDLQQLWKQNPKDIHTLAQ LISAYSLVDPEKAKALSKHLPSSDSMSLKV DVEALENSAGATYIRKKGGKVTGD SQPKEQGQGD LKKKKKKKKGKL PKNYDPKVTPDPERWLP MRERSYYRGRKK GKKKDQIGKGTQGATAGASSELDASKTVSSPPTSPRPGSAATVSASTSNIIPPR HQPAGAPATKKKQQQKKKKGGKGGW
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	SRP72
Protein Names	Recommended name: Signal recognition particle 72 kDa protein Short name= SRP72
Expression Region	2-671aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged



Mol. Weight	101.1 kDa
Protein Length	Full Length of Mature Protein

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

The process for recombinant Human SRP72 protein synthesis in e.coli cells includes recombining the DNA fragment that encodes the Human SRP72 protein (2-671aa) into expression vector, transforming the recombinant vector into e.coli cells, screening positive cells, cultivating and inducing cells, conducting cell lysis, and performing expression analysis. The resulting recombinant Human SRP72 protein is purified from the cell lysate using affinity purification, with its purity exceeding 85% as determined by SDS-PAGE.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.