




Recombinant Human Stress-induced-phosphoprotein 1 (STIP1)

Product Code	CSB-EP022831HU
Relevance	Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB).
Abbreviation	Recombinant Human STIP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31948
Alias	Hsc70/Hsp90-organizing protein ;HopRenal carcinoma antigen NY-REN-11Transformation-sensitive protein IEF SSP 3521
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MEQVNELKEKGNKALSVGNIDDALQCYSEAIKLDPHNHVLVSNRSAAYAKKGD YQKAYEDGCKTVDLKPDWVGKGYSRKAAALEFLNRFEEAKRTYEEGLKHEANN PQLKEGLQNMEARLAERKFMNPFNMPNLYQKLESDPRTRTLSDPTYRELIEQ LRNKPSDLGTLQDPRIMTTLSVLLGVDLGSMDEEEEIATPPPPPPPKKETKPE PMEEDLPENKKQALKEKELGNDAYKKKDFDTALKHYDKAKELDPTNMTYITNQ AAVYFEKGDYINKCRELCEKAIEVGRENREDYRQIAKAYARIGNSYFKEEKYKD AIHFYNKSLAEHRTPDVLKKCQQAEEKILKEQERLAYINPDALAEKNKGNECFQ KGDYPQAMKHYTEAIKRNPDAKLYSNRAACYTKLLEFQLALKDCEEICQLEPT FIKGYTRKAAALEAMKDYTEKAMDVYQKALDLDSSCKEAADGYQRCMMAQYN RHDSPEDVKRRAMADPEVQQIMSDPAMRLILEQMMDKDPQALSEHLKNPVIAQ KIQKLMDVGIIAIR
Research Area	Signal Transduction
Source	E.coli
Target Names	STIP1
Expression Region	1-543aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	78.6kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.