




Recombinant Human Sulfotransferase family cytosolic 2B member 1 (SULT2B1)

Product Code	CSB-EP022943HU
Relevance	Sulfotransferase that utilizes 3'-phospho-5'-adenylyl sulfate (PAPS) as sulfonate donor to catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs and xenobiotic compounds. Sulfonation increases the water solubility of most compounds, and therefore their renal excretion, but it can also result in bioactivation to form active metabolites. Sulfates hydroxysteroids like DHEA. Isoform 1 preferentially sulfonates cholesterol, and isoform 2 avidly sulfonates pregnenolone but not cholesterol.
Abbreviation	Recombinant Human SULT2B1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O00204
Alias	Alcohol sulfotransferase Hydroxysteroid sulfotransferase 2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MDGPAEPQIPGLWDTYEDDISEISQKLPGEYFRYKGVFPFVGLYSLESISLAEN TQDVRDDDIFIITYPKSGTTWMIEIICLILKEGDPSPWIRSVPIWERAPWCETIVGA FSLPDQYSPRLMSSHLPIQIFTKAFFSSKAKVIYMGRNPRDVVVSLYHYSKIAG QLKDPGTPDQFLRDLKGEVQFGSWFDHIKGWLRMGKNDNFLFITYEELQQD LQGSVERICGFLGRPLGKEALGSVVAHSTFSAMKANTMSNYTLLPPSLLDHRR GAFLRKGVCGDWKNHFTVAQSEAFDRAYRKQMRGMPTFPWDEDPPEEDGSP DPEPSPEPEPKPSLEPNTSLEREPRPNSSPSPSPGQASETPHPRPS
Research Area	Cardiovascular
Source	E.coli
Target Names	SULT2B1
Protein Names	Recommended name: Sulfotransferase family cytosolic 2B member 1 Short name= ST2B1 Short name= Sulfotransferase 2B1 EC= 2.8.2.2 Alternative name(s): Alcohol sulfotransferase Hydroxysteroid sulfotransferase 2
Expression Region	1-365aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	68.3kDa


Protein Length

Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.