




Recombinant Human Tescalcin (TESC)

Product Code	CSB-EP023393HU
Relevance	Functions as an integral cofactor in cell pH regulation by controlling plasma mbrane-type Na ⁺ /H ⁺ exchange activity. Promotes the maturation, transport, cell surface stability and exchange activity of SLC9A1/NHE1 at the plasma mbrane. Promotes the induction of hatopoietic st cell differentiation toward megakaryocytic lineage. Essential for the coupling of ERK cascade activation with the expression of ETS family genes in megakaryocytic differentiation. Also involved in granulocytic differentiation in a ERK-dependent manner. Inhibits the phosphatase activity of calcineurin.
Abbreviation	Recombinant Human TESC protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BS2
Alias	Tescalcin ;TSC
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	GAAHSASEEVRELEGKTFSSDQIEQLHRRFKQLSGDQPTIRKENFNNVPDLE LNPIRSKIVRAFFDNRNLRKGPSGLADEINFEDFLTIMSYFRPIDTTMDEEQVEL SRKEKLRFLFHMYSDDSDGRITLEEYRNVEELLSGNPHIEKESARSADGAM MEAASVCMGQMEPDQVYEGITFEDFLKIWQGIDIETKMHVRFNLNMETMALCH
Research Area	Signal Transduction
Source	E.coli
Target Names	TESC
Expression Region	2-214aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	40.6kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.