




Recombinant Human Cell surface hyaluronidase (TMEM2), partial

Product Code	CSB-EP023791HU
Abbreviation	Recombinant Human TMEM2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHN6
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SSKYAPDENCPCDQNPRLRNWDPGQDSAKQVVIKEGDMRLRLTSDATVHSIVIQ DGGLLVFGDNKDGSRNITLRTHYILIQDGGALHIGAEEKCRYKSKATITLYGKSDE GESMPTFGKKFIGVEAGGTLELHGARKASWTLLARTLNSS
Research Area	Cell Biology
Source	E.coli
Target Names	TMEM2
Protein Names	Recommended name: Transmembrane protein 2
Expression Region	104-250aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	21.5 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.