




# Recombinant Human Tubulin polymerization-promoting protein (TPPP)

<b>Product Code</b>	CSB-EP024115HU
<b>Relevance</b>	May play a role in the polymerization of tubulin into microtubules, microtubule bundling and the stabilization of existing microtubules, thus maintaining the integrity of the microtubule network. May play a role in mitotic spindle assembly and nuclear envelope breakdown.
<b>Abbreviation</b>	Recombinant Human TPPP protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O94811
<b>Alias</b>	25 kDa brain-specific protein TPPP/p25 p24 p25-alpha
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MADKAKPAKAANRTPPKSPGDPSKDRAAKRLSLESEGAGEGAAASPELSALE EAFRRFAVHGDA RATGREMHGKNWSKLCKDCQVIDGRNVTVDVDIVFSKIK GKSCRTITFEQFQEAL EELAKKRFKDKSSEEAVRE VHR LIEGKAPIISGVTKAIS SPTVSRLTDTTKFTGSHKERFDPSGKGKGKAGRVDLVDESGYVSGYKHAGTY DQKVQGGK
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	TPPP
<b>Protein Names</b>	Recommended name: Tubulin polymerization-promoting protein Short name= TPPP Alternative name(s): 25 kDa brain-specific protein TPPP/p25 p24 p25-alpha
<b>Expression Region</b>	1-219aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	50.7kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.