




Recombinant Human Splicing factor U2AF 35 kDa subunit (U2AF1)

Product Code	CSB-EP025405HU
Relevance	Plays a critical role in both constitutive and enhancer-dependent splicing by mediating protein-protein interactions and protein-RNA interactions required for accurate 3'-splice site selection. Recruits U2 snRNP to the branch point. Directly mediates interactions between U2AF2 and proteins bound to the enhancers and thus may function as a bridge between U2AF2 and the enhancer complex to recruit it to the adjacent intron.
Abbreviation	Recombinant Human U2AF1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q01081
Alias	U2 auxiliary factor 35 kDa subunit U2 small nuclear RNA auxiliary factor 1 U2 snRNP auxiliary factor small subunit
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAEYLASIFGTEKDKVNC SFYFKIGACRHGDRCSRLHNKPTFSQTIALLN IYRN PQNSSQSADGLRCAVSDVEMQEHYDEFFEEVFTEMEEEKYGEVEEMNVCDNL GDHLVGNVYVKFRREEDA EKAVIDLNNRW FNGQPIHAELSPVTD FREACCRQ YEMGECTRGGFCNFMHLKPISREL RRELYGRRRK KHRSRSR SRERRSRSRD RGRGGGGGGGGGGGGGRERDRRRSRDRERSGRF
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	U2AF1
Expression Region	1-240aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	54.9kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Description

Recombination of a plasmid encoding the Human U2AF1 protein (1-240aa) is the first step during the production the recombinant Human U2AF1 protein. The constructed plasmid is introduced into e.coli cells. e.coli cells that can survive in the presence of a specific antibiotic are selected to be cultured for the induction of protein expression. The protein is equipped with a N-terminal GST tag. After expression, affinity purification is used to isolate and purify the recombinant Human U2AF1 protein from the cell lysate. Denaturing SDS-PAGE is then applied to resolve the resulting recombinant Human U2AF1 protein. Its purity exceeds 90%.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.