




Recombinant Pig Vascular endothelial growth factor A (VEGFA)

Product Code	CSB-EP025833PI
Relevance	Growth factor active in angiogenesis, vasculogenesis and endothelial cell growth. Induces endothelial cell proliferation, promotes cell migration, inhibits apoptosis and induces permeabilization of blood vessels. Binds to the FLT1/VEGFR1 and KDR/VEGFR2 receptors, heparan sulfate and heparin .
Abbreviation	Recombinant Pig VEGFA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49151
Product Type	Recombinant Proteins
Immunogen Species	Sus scrofa (Pig)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	APMAEGDQKPHEVVKFMDVYQRSYCRPIETLVDIFQEYPDEIEYIFKPSCVPLM RCGGCCNDEGLECVPTTEEFNITMQIMRIKPHQGQHIGEMSFLQHNKCECRPK KDRARQENPCGPCSERRKHLFVQDPQTCKCCKNTDSRCKARQLELNERTC RCDKPRR
Research Area	Others
Source	E.coli
Target Names	VEGFA
Expression Region	27-190aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	23.2kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.