




Recombinant Rabbit Vascular endothelial growth factor A (VEGFA)

Product Code	CSB-EP025833RB
Abbreviation	Recombinant Rabbit VEGFA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	XP_002714743.1
Product Type	Recombinant Proteins
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	YLWLLHAPLLHVSLDSLVAAFSVVFEVLPSSLDIPGSKKEEQASTQLGGVEEEH EASGVVTRQLAFVVDITNPTLEKETKTITRIHKAPKFPSHFGFWKRAEEERGG KSSGEKSRKTDGVGERAQAGEQREGPGQRAAALTDRQTDAPSPSAHLLPG RRPTVDAAASGGQQPEPAPEGGVEGVGARGIALKLFVQLLGCSRSGGVVRA GGAEPSCAARSVSSGREPPPPPPPEEEGEEEGEKEEEERGPRWRLGAGEPGS WTGEAAVCADSAPAARAPQALARASAPGARGARGGAEESGPSRSPSRRGSA SRAGPGRASETMNFLLSVHWSLALLLYLHHAKWSQAAPMAEEGDNKPHEV VKFMEVYRRSYCQPIETLVDIFQEYPDEIEYIFKPSCVPLVRCGGCCNDESLEC VPTEEFNVMTQMIRIKPHQGQHIGEMSFLQHKNKCECRCDKPRR
Research Area	Others
Source	E.coli
Target Names	VEGFA
Expression Region	51-511aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	65.7kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.