




Recombinant Bacillus subtilis Uncharacterized protein yuaB (yuaB)

Product Code	CSB-EP300296BRJ
Relevance	Required for complex colony architecture.
Abbreviation	Recombinant Bacillus subtilis bsIA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P71014
Product Type	Recombinant Proteins
Immunogen Species	Bacillus subtilis (strain 168)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MKRKLLSSLAISLSLGLLVSAPTASFAAESTSTKAHTESTMRTQSTASLFATIT GASKTEWSFSDIELTYRPNTLLSLGVMEFTLP SGFTANTKDTLNGNALRTTQIL NNGKTVRVPLALDLLGAGEFKLKLNNKTLPAAGTYTFRAENKSLSIGNKFYAEA SIDVAKRSTPPTQPCGCN
Research Area	Microbiology
Source	E.coli
Target Names	yuaB
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	35.3kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.