




Recombinant Escherichia coli Oxidoreductase YdhF (ydhF)

Product Code	CSB-EP300544ENV
Relevance	May function as oxidoreductase.
Abbreviation	Recombinant E.coli ydhF protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P76187
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MVQRITIA PQGPEFSRFVMGYWRLMDWNMSARQLVSFIEEHLDLGVTTVDHA DIYGGYQCEAAFGEALKLAPHLRERMEIVSKCGIATTAREENVIGHYITDRDHII KSAEQSLINLATDHLDLLLIHRPDPLMDADEVADAFKHLHQSGKVRHFGVSNFT PAQFALLQSRLPFTLATNQVEISPVHQPLLLDGTLDQLQQLRVRPMAWSCCLGG GRLFNDYFQPLRDELAVVAEELNAGSIEQVVYAWVLRLPSQPLPIIGSGKIER VRAAVEAETLKMTRQQWFRIRKAALGYDVP
Research Area	Others
Source	E.coli
Target Names	ydhF
Expression Region	1-298aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	37.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.