




Recombinant Vaccinia virus Protein A33 (A33R), partial

Product Code	CSB-EP300755VAA1
Abbreviation	Recombinant Vaccinia virus A33R protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P68616
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	VRLNQCMSANEAAITDAAVAVAAASSTHRKVASSTTQYDHKESCNGLYYQGS CYILHSDYQLFSDAKANCTAESSTLPNKSVDLITWLIDYVEDTWGSDGNPITKT TSDYQSDVSQEVRYFCVKTMN
Research Area	Others
Source	E.coli
Target Names	A33R
Protein Names	Recommended name: Protein A33
Expression Region	57-185aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	21.7 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.