




Recombinant Acetoanaerobium sticklandii Ferredoxin (CLOST_2292)

Product Code	CSB-EP301302DUS
Abbreviation	Recombinant Acetoanaerobium sticklandii CLOST_2292 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P80168
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Acetoanaerobium sticklandii (strain ATCC 12662 / DSM 519 / JCM 1433 / NCIMB 10654) (Clostridium sticklandii)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	AYVINDSCISCGACEPECPVNAITAGDDKYVIDAATCIDCGACAGVCPVDAPQP E
Research Area	Others
Source	E.coli
Target Names	CLOST_2292
Protein Names	Recommended name: Ferredoxin
Expression Region	2-56aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	13.0 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.