




Recombinant Enterobacteria phage f1 Replication-associated protein G2P (II)

Product Code	CSB-EP304002ECHb1
Abbreviation	Recombinant Enterobacteria phage f1 Replication-associated protein G2P
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P69546
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Enterobacteria phage f1 Replication-associated protein G2P(II)
Immunogen Species	Enterobacteria phage f1 (Bacteriophage f1)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MIDMLVLRLPFIDSLVCSRLSGNDLIAFVDLSKIATLSGMNLSARTVEYHIDGDLT VSGLSHPFESLPHTHYSGIAFKIYEGSKNFYPCVEIKASPAKVLQGHNVFGTTDL ALCSEALLLNFANSLPCLYDLLDVNATTISRIDATFSARAPNENIAKQVIDHLRNV SNGQTKSTRSQNWESTVTWNETSRRHRTL VAYLKHVELQHQQQLSSKPSAKM TSYQKEQLKVLSNPDLLEFASGLVRFEARIETRYLKSFGPLPLNLFDAIRFASDYN SQGKDLIFDLWSFSFSELFKA FEGDSMNIYDDSAVLDAIQSKHFTITPSGKTSFA KASRYFGFYRRLVNEGYDSVALTMPRNSFWRYVSALVECGIPKSQLMNLSTC NNVPLVRFINVDFSSQRPDWAYNEPVLKIA
Research Area	Others
Source	E.coli
Target Names	II
Expression Region	1-410aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	53.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.