




Recombinant *Gossypium hirsutum* zinc finger protein 8-like (LOC107899218)

Product Code	CSB-EP3066GHB
Abbreviation	Recombinant <i>Gossypium hirsutum</i> LOC107899218 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A1U8IQM3
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MEKNERETHDFMNVESFSQLPFIRPAPSKEKGIRLFGKEFGGVDPATPSNESD SAENNEDTTKENENNGDNSRRFECHYCCRNFPSTSQALGGHQNAHKRERQHA KRAHLQSAMVHTSLSDAHYGLVNYRLGSAPTPPITYPSWNSSFTRSTSRFYG NHTSFSHHPINGSPLGLWRIPSTLQNNSSNFNPDRSSSSSSSHPLPLFAGDE LKPPSQVVAGGGGSSSQSRYVYESKPRLQDHVSLDLHL
Research Area	Others
Source	E.coli
Target Names	LOC107899218
Expression Region	1-248aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	31.7 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.