




Recombinant *Saccharomyces cerevisiae* Ubiquitin-like-specific protease 1 (ULP1), partial

Product Code	CSB-EP311619SVG
Abbreviation	Recombinant <i>Saccharomyces cerevisiae</i> ULP1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q02724
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	LVPRGSHMASLVPENLKDDDDQVQKALASRENTQLMNRDNIEITVRDFKTLAP RRWLNDTIIIEFFMKYIEKSTPNTVAFNSFFYTNLSERGYQGVRWRMKRKKTKI DKLDKIFTPINLNQSHWALGIIDLKKKTIGYVDSLNGPNAMSFALITDLQKYVM EESKHTIGEDFDLIHLDCPQQPNGYDCGIYVCMNTLYGSADAPLDFDYKDAIR MRRFIAHLILTDAK
Research Area	Others
Source	E.coli
Target Names	ULP1
Protein Names	Recommended name: Ubiquitin-like-specific protease 1 EC= 3.4.22.68
Expression Region	403-621aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	27.9 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.