




Recombinant Salmonella typhimurium Invasion protein iagB (iagB)

Product Code	CSB-EP315647SXB
Abbreviation	Recombinant Salmonella typhimurium iagB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0CL15
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	DCWLQAEKMFNIESELLYAIAQQESAMKPGAIGHNRDGSTDLGLMQINSFHMK RLKKMGISEKQLLQDPCISVIVGASILSDMMKIYGYSWEAVGAYNAGTSPKRSD IRKRYAKKIWENYRKLKGMSAEKKNKRLSIAANK
Research Area	Others
Source	E.coli
Target Names	iagB
Expression Region	20-160aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	32.0kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.